

2,3,3',5,5'-Pentachlorobiphenyl 100 µg/mL in Isooctane

Art. ID	CIL-PCB-111-CS
Unit	1,2 mL
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,3,3',5,5'-Pentachloro biphenyl (PCB 111)	[39635-32-0]	100	µg/mL		