

## **2',3,4,4',5-Pentachlorobiphenyl 100 µg/mL in Isooctane**

Art. ID	CIL-PCB-123-CS
Unit	1,2 mL
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

### Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2',3,4,4',5-Pentachloro biphenyl (PCB 123)	[65510-44-3]	100	µg/mL		