

## **4,4'-Dichlorobiphenyl 100 µg/mL in Isooctane**

|                 |  |
|-----------------|--|
| Art. ID         | CIL-PCB-15-CS  |
| Unit            | 1,2 mL   |
| Deliverydetails | ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2 |

### Description

Chemical purity: 98%

| Text/Information | Analyte/Parameter                  | CAS number  | Concentration/Value | Unit  | Method | Source |
|------------------|------------------------------------|-------------|---------------------|-------|--------|--------|
|                  | 4,4'-Dichlorobiphenyl (<br>PCB 15) | [2050-68-2] | 100                 | µg/mL |        |        |