

### **4,4'-Dichlorobiphenyl 100 µg/mL in Isooctane**

Art. ID	CIL-PCB-15-CS
Unit	1,2 mL
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

#### Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4,4'-Dichlorobiphenyl ( PCB 15)	[2050-68-2]	100	µg/mL		