

2,3,3',4,4',5'-Hexachlorobiphenyl 35 µg/mL in Isooctane

Art. ID CIL-PCB-157
Unit 1 mL
Deliverydetails De Minimis UN1262 II E2

Description

Chemical purity: 98%

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---|--------------|---------------------|-------|--------|--------|
| | 2,3,3',4,4',5'-Hexachlorobiphenyl (PCB 157) | [69782-90-7] | 35 | µg/mL | | |