

## **2,2',3,3',4,4',5-Heptachlorobiphenyl 100 µg/mL in Isooctane**

Art. ID CIL-PCB-170-CS  
Unit 1,2 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

### Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,2',3,3',4,4',5-Heptachlorobiphenyl (PCB 170)	[35065-30-6]	100	µg/mL		