

2,2',3,3',4',5,6-Heptachlorobiphenyl 100 µg/mL in Isooctane

Art. ID CIL-PCB-177-CS
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,2',3,3',4',5,6-Heptachlorobiphenyl (PCB 177)	[52663-70-4]	100	µg/mL		