

2,3,4-Trichlorobiphenyl 35 µg/mL in Isooctane

| | |
|-----------------|-------------------------|
| Art. ID | CIL-PCB-21 |
| Unit | 1 mL |
| Deliverydetails | De Minimis UN1262 II E2 |

Description

Chemical purity: 98%

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------------------------|--------------|---------------------|-------|--------|--------|
| | 2,3,4-Trichlorobiphenyl (PCB 21) | [55702-46-0] | 35 | µg/mL | | |