

## **2,3,4'-Trichlorobiphenyl 35 µg/mL in Isooctane**

|                 |                         |
|-----------------|-------------------------|
| Art. ID         | CIL-PCB-22              |
| Unit            | 1 mL                    |
| Deliverydetails | De Minimis UN1262 II E2 |

### Description

Chemical purity: 98%

| Text/Information | Analyte/Parameter                      | CAS number   | Concentration/Value | Unit  | Method | Source |
|------------------|--|--------------|---------------------|-------|--------|--------|
|                  | 2,3,4'-Trichlorobiphenyl<br>I (PCB 22) | [38444-85-8] | 35                  | µg/mL |        |        |