

### **3,4,4'-Trichlorobiphenyl 100 µg/mL in Isooctane**

Art. ID	CIL-PCB-37-CS
Unit	1,2 mL
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

#### Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3,4,4'-Trichlorobiphenyl I (PCB 37)	[38444-90-5]	100	µg/mL		