

2,2',4,5'-Tetrachlorobiphenyl 100 µg/mL in Isooctane

Art. ID CIL-PCB-49-CS
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

Description

Chemical purity: 98%

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--|--------------|---------------------|-------|--------|--------|
| | 2,2',4,5'-Tetrachlorobiphenyl (PCB 49) | [41464-40-8] | 100 | µg/mL | | |