

Phosmet (unlabelled) 100 µg/mL in Acetonitrile

Art. ID CIL-ULM-8454-1.2
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Description

Chemical purity: 98%

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|-------|--------|--------|
| | Phosmet | [732-11-6] | 100 | µg/mL | | |