

CONOSTAN - CONOSTAN Custom Blend AM-special in 75 blank oil: AM Special 300 ppm

Art. ID CON-150-250-008
Unit 100 g
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Barium (Ba)	[7440-39-3]	0,03	%		
	Barium (Ba)	[7440-39-3]	300	µg/g		
	Calcium (Ca)	[7440-70-2]	0,03	%		
	Calcium (Ca)	[7440-70-2]	300	µg/g		
	Magnesium (Mg)	[7439-95-4]	0,03	%		
	Magnesium (Mg)	[7439-95-4]	300	µg/g		
	Phosphorus (P)	[7723-14-0]	0,03	%		
	Phosphorus (P)	[7723-14-0]	300	µg/g		
	Zinc (Zn)	[7440-66-6]	0,03	%		
	Zinc (Zn)	[7440-66-6]	300	µg/g		