

Gymnodimine (GYM) standard solution in Methanol

Art. ID	NRC-CRM-GYM-b
Unit	0,5 mL
Deliverydetails	De Minimis UN1230 (6.1) II E2

Description

NRC-CRM-GYM-b is a certified calibration solution for the determination of gymnodimine (GYM) by chemical analysis methods. Each ampoule contains approx. 0.5 mL of a solution of GYM in methanol at a concentration of 4.93 µM. GYM is a member of the cyclic imine group of toxins.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Gymnodimine (GYM)		4,93	µM		