

**Arsenic Standard for ICP/MS - Claritas PPT Grade - Arsenic (As) 1 mg/L in 2% HNO<sub>3</sub>**

Art. ID LM24-SP-CLAS2-1BY  
Unit 125 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Arsenic (As)      | [7440-38-2] | 1                   | mg/L |        |        |