

Aluminum Organometallic Oil Standard for Wear Metal Analysis - Aluminium (Al) 5000 µg/g in Base Oil 75

Art. ID LM24-SP-ORG-AL8-4Z
Unit 50 g
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Aluminium (Al) | [7429-90-5] | 5000 | µg/g | | |