

Beryllium Organometallic Oil Standard for Wear Metal Analysis - Beryllium (Be) 1000 µg/g in Base Oil 75

Art. ID LM24-SP-ORG-BE8-2Z
Unit 50 g
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Beryllium (Be) | [7440-41-7] | 1000 | µg/g | | |