

Copper Organometallic Oil Standard for Wear Metal Analysis - Copper (Cu) 1000 µg/g in Base Oil 75

Art. ID LM24-SP-ORG-CU8-2Z
Unit 50 g
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Copper (Cu) | [7440-50-8] | 1000 | µg/g | | |