

Tin Organometallic Oil Standard for Wear Metal Analysis - T in (Sn) 5000 µg/g in Base Oil 75

Art. ID LM24-SP-ORG-SN8-4Z
Unit 50 g
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Tin (Sn) | [7440-31-5] | 5000 | µg/g | | |