

**Assurance Grade Arsenic Standard for AA and ICP - Arsenic (As) 10000 mg/L in 5% HNO<sub>3</sub>**

Art. ID LM24-SP-PLAS2-3Y  
Unit 125 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	10000	mg/L		