

Assurance Grade Copper Standard for AA and ICP - Copper (Cu) 1000 mg/L in 2% HNO₃

Art. ID LM24-SP-PLCU2-2T
Unit 250 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Copper (Cu) | [7440-50-8] | 1000 | mg/L | | |