

Assurance Grade Antimony Standard for AA and ICP - Antimony (Sb) 10000 mg/L in Water/6% Tartaric acid/1% HNO₃

Art. ID LM24-SP-PLSB7-3X
Unit 500 mL
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Antimony (Sb) | [7440-36-0] | 10000 | mg/L | | |