

Assurance Grade Selenium Standard for AA and ICP - Selenium (Se) 10000 mg/L in 5% HNO₃

Art. ID LM24-SP-PLSE2-3Y
Unit 125 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Selenium (Se) | [7782-49-2] | 10000 | mg/L | | |