

Assurance Grade Tin Standard for AA and ICP - T in (Sn) 1000 mg/L in 20% HCl

Art. ID LM24-SP-PLSN5-2T
Unit 250 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1789 II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Tin (Sn) | [7440-31-5] | 1000 | mg/L | | |