

Assurance Grade Tungsten Standard for AA and ICP - Tungsten (W) 10000 mg/L in 2% HNO₃/5% HF

Art. ID LM24-SP-PLW2-3X
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Tungsten (W) | [7440-33-7] | 10000 | mg/L | | |