

**Assurance Grade Tungsten Standard for AA and ICP - Tungsten (W) 1000 mg/L in Water**

Art. ID LM24-SP-PLW9-2Y  
Unit 125 mL  
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Tungsten (W)      | [7440-33-7] | 1000                | mg/L |        |        |