

**Assurance Grade Tungsten Standard for AA and ICP - Tungsten (W) 10000 mg/L in Water**

Art. ID LM24-SP-PLW9-3X  
Unit 500 mL  
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Tungsten (W)      | [7440-33-7] | 10000               | mg/L |        |        |