

2,2',5,5'-Tetrachlorobiphenyl (PCB 52) 50 µg/mL in Acetonitrile

Art. ID LM24-SP-QUECR-PCB-52

Unit 5 mL

Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,2',5,5'-Tetrachlorobi phenyl (PCB 52)	[35693-99-3]	50	µg/mL		