

Assurance Grade instrument Check (Lab Performance) Standard 4 in 2% HNO₃

Art. ID LM24-SP-CALMIX4-100
Unit 125 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Barium (Ba)	[7440-39-3]	1	mg/L		
	Aluminium (Al)	[7429-90-5]	10	mg/L		
	Arsenic (As)	[7440-38-2]	10	mg/L		
	Copper (Cu)	[7440-50-8]	10	mg/L		
	Manganese (Mn)	[7439-96-5]	10	mg/L		
	Sodium (Na)	[7440-23-5]	10	mg/L		
	Nickel (Ni)	[7440-02-0]	10	mg/L		
	Phosphorus (P)	[7723-14-0]	10	mg/L		
	Lead (Pb)	[7439-92-1]	10	mg/L		
	Scandium (Sc)	[7440-20-2]	10	mg/L		
	Zinc (Zn)	[7440-66-6]	10	mg/L		
	Potassium (K)	[7440-09-7]	50	mg/L		