

ICP-MS interference A2 Standard in 5% HNO₃/tr HF-125 ML in 5% HNO₃/tr. HF

Art. ID LM24-SP-CL-INT-A3

Unit 125 mL

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Molybdenum (Mo)	[7439-98-7]	20	mg/L		
	Titanium (Ti)	[7440-32-6]	20	mg/L		
	Aluminium (Al)	[7429-90-5]	1000	mg/L		
	Potassium (K)	[7440-09-7]	1000	mg/L		
	Magnesium (Mg)	[7439-95-4]	1000	mg/L		
	Phosphorus (P)	[7723-14-0]	1000	mg/L		
	Sulfur (S)	[7704-34-9]	1000	mg/L		
	Carbon (C)	[7440-44-0]	2000	mg/L		
	Iron (Fe)	[7439-89-6]	2500	mg/L		
	Sodium (Na)	[7440-23-5]	2500	mg/L		
	Calcium (Ca)	[7440-70-2]	3000	mg/L		
	Chlorine (Cl)	[7782-50-5]	20000	mg/L		