

### **Assurance Grade instrument Calibration Standard 3 in 5% HNO<sub>3</sub>**

Art. ID LM24-SP-ICAL-3  
Unit 125 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	50	mg/L		
	Copper (Cu)	[7440-50-8]	250	mg/L		
	Cobalt (Co)	[7440-48-4]	500	mg/L		
	Vanadium (V)	[7440-62-2]	500	mg/L		
	Iron (Fe)	[7439-89-6]	1000	mg/L		
	Aluminium (Al)	[7429-90-5]	2000	mg/L		
	Barium (Ba)	[7440-39-3]	2000	mg/L		