

### **Assurance Grade instrument Calibration Standard 4A in 5% HNO<sub>3</sub>**

Art. ID LM24-SP-ICAL-4A

Unit 125 mL

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Lead (Pb)	[7439-92-1]	30	mg/L		
	Cadmium (Cd)	[7440-43-9]	50	mg/L		
	Selenium (Se)	[7782-49-2]	50	mg/L		
	Arsenic (As)	[7440-38-2]	100	mg/L		
	Thallium (Tl)	[7440-28-0]	100	mg/L		