

instrument Check Standard 3 for Ion Chromatography in 2% HNO3

Art. ID LM24-SP-ICMIX3-100
Unit 125 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Li+		50	mg/L		
	K+		200	mg/L		
	Mg+2		200	mg/L		
	Na+		200	mg/L		
	Ammonium (NH ₄ ⁺)	[14798-03-9]	400	mg/L		
	Ca+2		1000	mg/L		