

instrument Check Standard 4 for Ion Chromatography in 0.5% HNO₃

Art. ID LM24-SP-ICMIX4-100
Unit 125 mL
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|------------------------------------------|--------------|---------------------|------|--------|--------|
| | Li+ | | 10 | mg/L | | |
| | Na+ | | 50 | mg/L | | |
| | K+ | | 100 | mg/L | | |
| | Ammonium (NH ₄ ⁺) | [14798-03-9] | 100 | mg/L | | |