

### **Assurance Grade Mixed Multi-Element Calibration Standard 1 in 2% HNO<sub>3</sub>**

Art. ID LM24-SP-MIXSTD1-500  
Unit 500 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	50	mg/L		
	Manganese (Mn)	[7439-96-5]	100	mg/L		
	Cadmium (Cd)	[7440-43-9]	150	mg/L		
	Zinc (Zn)	[7440-66-6]	150	mg/L		
	Selenium (Se)	[7782-49-2]	200	mg/L		
	Lead (Pb)	[7439-92-1]	500	mg/L		