

Assurance Grade Mixed Multi-Element Calibration Standard 1A in 5% HNO₃/Tr. Tartaric acid/Tr. HF

Art. ID LM24-SP-MIXSTD1A-100
Unit 125 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Silver (Ag) | [7440-22-4] | 5 | mg/L | | |
| | Barium (Ba) | [7440-39-3] | 10 | mg/L | | |
| | Boron (B) | [7440-42-8] | 20 | mg/L | | |
| | Cadmium (Cd) | [7440-43-9] | 20 | mg/L | | |
| | Copper (Cu) | [7440-50-8] | 20 | mg/L | | |
| | Manganese (Mn) | [7439-96-5] | 20 | mg/L | | |
| | Antimony (Sb) | [7440-36-0] | 50 | mg/L | | |
| | Selenium (Se) | [7782-49-2] | 50 | mg/L | | |
| | Arsenic (As) | [7440-38-2] | 100 | mg/L | | |
| | Calcium (Ca) | [7440-70-2] | 100 | mg/L | | |