

AA/ICP Mixed Multi-Element Calibration Standard 1C, 500 mL in 5% HNO₃/Tr. Tartaric acid/Tr. HF

Art. ID LM24-SP-MIXSTD1C-500
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	5	mg/L		
	Boron (B)	[7440-42-8]	10	mg/L		
	Barium (Ba)	[7440-39-3]	10	mg/L		
	Cadmium (Cd)	[7440-43-9]	20	mg/L		
	Copper (Cu)	[7440-50-8]	20	mg/L		
	Manganese (Mn)	[7439-96-5]	20	mg/L		
	Antimony (Sb)	[7440-36-0]	50	mg/L		
	Selenium (Se)	[7782-49-2]	50	mg/L		
	Arsenic (As)	[7440-38-2]	100	mg/L		
	Calcium (Ca)	[7440-70-2]	100	mg/L		