

Assurance Grade Mixed Multi-Element Calibration Standard 2 in 5% HNO₃

Art. ID LM24-SP-MIXSTD2-500
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Barium (Ba)	[7440-39-3]	100	mg/L		
	Cobalt (Co)	[7440-48-4]	100	mg/L		
	Copper (Cu)	[7440-50-8]	100	mg/L		
	Vanadium (V)	[7440-62-2]	100	mg/L		
	Iron (Fe)	[7439-89-6]	10000	mg/L		