

### **Assurance Grade Mixed Multi-Element Calibration Standard 2A in 5% HNO<sub>3</sub>**

Art. ID LM24-SP-MIXSTD2A-500  
Unit 500 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Strontium (Sr)	[7440-24-6]	10	mg/L		
	Lithium (Li)	[7439-93-2]	50	mg/L		
	Molybdenum (Mo)	[7439-98-7]	100	mg/L		
	Sodium (Na)	[7440-23-5]	100	mg/L		
	Potassium (K)	[7440-09-7]	200	mg/L		