

Assurance Grade Mixed Multi-Element Calibration Standard 3 in 2% HNO₃/Tr. HF

Art. ID LM24-SP-MIXSTD3-100
Unit 125 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Molybdenum (Mo)	[7439-98-7]	100	mg/L		
	Silicon (Si)	[7440-21-3]	100	mg/L		
	Arsenic (As)	[7440-38-2]	500	mg/L		