

Assurance Grade Mixed Multi-Element Calibration Standard 3 in 2% HNO₃/Tr. HF

Art. ID LM24-SP-MIXSTD3-500
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Molybdenum (Mo) | [7439-98-7] | 100 | mg/L | | |
| | Silicon (Si) | [7440-21-3] | 100 | mg/L | | |
| | Arsenic (As) | [7440-38-2] | 500 | mg/L | | |