

Assurance Grade Mixed Multi-Element Calibration Standard 4 in 5% HNO₃

Art. ID LM24-SP-MIXSTD4-500

Unit 500 mL

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Chromium (Cr)	[7440-47-3]	20	mg/L		
	Nickel (Ni)	[7440-02-0]	20	mg/L		
	Aluminium (Al)	[7429-90-5]	200	mg/L		
	Sodium (Na)	[7440-23-5]	200	mg/L		
	Potassium (K)	[7440-09-7]	400	mg/L		
	Calcium (Ca)	[7440-70-2]	1000	mg/L		