

Assurance Grade Mixed Multi-Element Calibration Standard 5 in 5% HNO₃/tr. Tartaric acid/tr. HF

Art. ID LM24-SP-MIXSTD5-100
Unit 125 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	50	mg/L		
	Boron (B)	[7440-42-8]	100	mg/L		
	Antimony (Sb)	[7440-36-0]	200	mg/L		
	Thallium (Tl)	[7440-28-0]	200	mg/L		
	Magnesium (Mg)	[7439-95-4]	1000	mg/L		