

Assurance Grade Mixed Multi-Element Calibration Standard 5A in 5% HNO₃

Art. ID LM24-SP-MIXSTD5A-500
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	10	mg/L		
	Nickel (Ni)	[7440-02-0]	20	mg/L		
	Thallium (Tl)	[7440-28-0]	50	mg/L		
	Iron (Fe)	[7439-89-6]	100	mg/L		
	Magnesium (Mg)	[7439-95-4]	100	mg/L		
	Lead (Pb)	[7439-92-1]	100	mg/L		