

**Assurance Grade Mixed Multi-Element Calibration Standard 4A, without Mercury in 5% HNO<sub>3</sub>/tr. HF**

Art. ID LM24-SP-MXSTD4A-100N  
Unit 125 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tin (Sn)	[7440-31-5]	40	mg/L		
	Chromium (Cr)	[7440-47-3]	50	mg/L		
	Zinc (Zn)	[7440-66-6]	50	mg/L		
	Aluminium (Al)	[7429-90-5]	100	mg/L		
	SiO <sub>2</sub>		100	mg/L		
	Titanium (Ti)	[7440-32-6]	100	mg/L		