

## **Assurance Grade Quality Control Standard 7A in 5% HNO<sub>3</sub>/Tr. HF**

Art. ID LM24-SP-QC-7A-500  
Unit 500 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	50	mg/L		
	Aluminium (Al)	[7429-90-5]	100	mg/L		
	Boron (B)	[7440-42-8]	100	mg/L		
	Barium (Ba)	[7440-39-3]	100	mg/L		
	Sodium (Na)	[7440-23-5]	100	mg/L		
	Silicon (Si)	[7440-21-3]	500	mg/L		
	Potassium (K)	[7440-09-7]	1000	mg/L		