

12 Element Organometallic Oil Standard for Wear Metal Analysis in Base Oil 75

Art. ID LM24-SP-S12-100Y
Unit 100 g
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	100	µg/g		
	Aluminium (Al)	[7429-90-5]	100	µg/g		
	Chromium (Cr)	[7440-47-3]	100	µg/g		
	Copper (Cu)	[7440-50-8]	100	µg/g		
	Iron (Fe)	[7439-89-6]	100	µg/g		
	Magnesium (Mg)	[7439-95-4]	100	µg/g		
	Sodium (Na)	[7440-23-5]	100	µg/g		
	Nickel (Ni)	[7440-02-0]	100	µg/g		
	Lead (Pb)	[7439-92-1]	100	µg/g		
	Silicon (Si)	[7440-21-3]	100	µg/g		
	Tin (Sn)	[7440-31-5]	100	µg/g		
	Titanium (Ti)	[7440-32-6]	100	µg/g		