

Semi-Volatile 4-Component GC/MS Tuning Standard in Methylene chloride (Dichloromethane)

Art. ID LM24-SP-CLPS-T4
Unit 1 mL
Deliverydetails De Minimis UN1593 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4,4'-Benzidine	[92-87-5]	2500	µg/mL		
	Bis(pentafluorophenyl)phenylphosphine	[5074-71-5]	2500	µg/mL		
	4,4'-DDT	[50-29-3]	2500	µg/mL		
	Pentachlorophenol	[87-86-5]	2500	µg/mL		