

QuEChERS GC/MS internal Standard Mix in Acetonitrile - 5 mL in Acetonitrile

Art. ID LM24-SP-QUECR-MS-IS

Unit 5 mL

Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Triphenylmethan	[519-73-3]	10	µg/mL		
	Triphenylphosphate	[115-86-6]	20	µg/mL		
	2,2',5,5'-Tetrachlorobi phenyl (PCB 52)	[35693-99-3]	50	µg/mL		
	2,2',5-Trichlorobipheny I (PCB 18)	[37680-65-2]	50	µg/mL		
	2,4,4'-Trichlorobipheny I (PCB 28)	[7012-37-5]	50	µg/mL		
	Tris(2-chloro-1-(chloro methyl)ethyl)phosphate	[13674-87-8]	50	µg/mL		