

QuEChERS QC Standard 1 in Acetonitrile - 5 mL in Acetonitrile

Art. ID LM24-SP-QUECR-MS-QC-1
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,2',3,4,4',5'-Hexachlorobiphenyl (PCB 138)	[35065-28-2]	50	µg/mL		
	2,2',4,4',5,5'-Hexachlorobiphenyl (PCB 153)	[35065-27-1]	50	µg/mL		